











Part number	Description	fits model	photo
A1470075	Soft seat overlay, hygiene recess, fits OceanVIP	OceanVIP	
A1470077	Arm locks	Ocean	
A10230	Replacement collection pan with lid	all Ocean models	
D10167	Safety bar, metal bar with foam grip surface, for all Ocean models except XL.	OceanVIP, Ocean, OceanSP	
A16511	Splash guard for Ocean models with hygiene soft seat	all Ocean models	
A16343	Key insert, fits hygiene soft seat opening to make solid seating surface	all Ocean models	
A1470086	Arm width extension	OceanVIP	
A16342	Soft seat toilet opening	all Ocean models	
A1470085	Reclining back support	OceanVIP	
A1470083	Amputee support, left	OceanVIP	
A1470027	Amputee support, right	Ocean VIP	

A1470084	Calf support, 2 pads, left	OceanVIP	
A1470828	Calf support, 2 pads, right	OceanVIP	
D100031	Action Gel seat	all Ocean models	
RP052611197	Side lateral, pair	all Ocean models	
RP055594	Arm troughs, pair, 12" long, recommend using arm width extension so arm will clear when flipped up.	all Ocean models	
A16521	Carry bag, Cordura Nylon, luggage wheels and shoulder strap.	Ocean, OceanSP and OceanXL	
D10175	Calf strap, 4" wide with hook and loop closure	all Ocean models	
A1470081	Padded chest strap, 4" wide with buckle closure	all Ocean models, included with EVIP	
A1470082	Padded lap belt, 4" wide with buckle closure	all Ocean models	
D10068	Safety harness, shoulder and lap belt combination	all Ocean models	

