























	Product number	Description	Diameter grab bar	Diameter suction pad	Grab length, (between pads)	Overall dimension	Capacity lbs	Product Weight lbs
	R1403999	Single hand grip, fixed length, 4" grab bar length	1.37"	4.69"	4"	5"W x 2.75"H x 12.75"L	150	1.5
	R1400221	Small bar telescopes from 8.25"-13.25"	1.37"	4.69"	8.5" – 13.4"	5"W x 2.75"H x 17" – 22"L	144	1.8
	R1400221S	Small bar telescopes from 8.25" - 13.25", with suction indicator on suction pad	1.37"	4.69"	8.5" – 13.4"	5"W x 2.75"W x 17" – 22"L	144	1.9
	R1400222	Medium bar telescopes from 13"-17.9"	1.37"	4.69"	13.2" – 18"	5"W x 2.75" H x 21.8" – 26.7"L	144	2.1
	R1400222S	Medium bar telescopes from 13"-17.9", with suction indicator on suction pad	1.37"	4.69"	13.2" – 18:	5"W x 2.75" H x 21.8" – 26.7"L	144	2.25
	R1400223	Large bar telescopes from 16.9"-21.8"	1.37"	4.69"	17" – 22"	5"W x 2.75" H x 25.7" – 30.6"L	144	2.4
	R1400223S	Large bar telescopes from 16.9"-21.8", with suction indicator on suction pad	1.37"	4.69"	17" – 22"	5"W x 2.75"H x 25.7" – 30.6"L	144	2.5
	R1400224	Xlarge bar telescopes from 26"-31"	1.37"	4.69"	25.9" – 30.8"	5"W x 2.75" H x 34.5" – 39.4"L	132	2.9
	R1400224S	Xlarge bar telescopes from 26"-31", with suction indicator on suction pad	1.37"	4.69"	25.9" – 30.8"	5"W x 2.75" H x 34.5" – 39.4"L	132	3
	R1400225	2 hand grip, fixed length 8.3"	1.37"	4.69"	8"	5"W x 2.75" H x 16.7"L	166	1.8
	R1400225S	2 hand grip 8.3" fixed length, with suction indicator on suction pad	1.37"	4.69"	8"	5"W x 2.75"H x 16.7"L	166	1.8
	R1400226	2 hand grip fixed length 13.8"	1.37"	4.69"	13.6"	5"W x 2.75"H x 22.2"L	154	2.1
	R1400226S	2 hand grip 13.8" with suction indicator button	1.37"	4.69"	13.6"	5" W x 2.75" H x 22.2"L	154	2.2
	R1400227	2 hand grip fixed length 25.6"	1.37"	4.69"	25.4"	5"W x 2.75"H x 33.9"L	154	2.8
	R1400227S	2 hand grip 25.6" with suction indicator button	1.37"	4.69"	25.4"	5" W X 2.75"H x 33.9"L	154	2.9

	R1400228	Kid Grip fixed length 6.7", thinner diameter grab bar	1"	4.69"	6.6"	5"W x 2.75"H x 15.2"L	100	1.6
	R1400228S	Kid Grip 6.7" with suction indicator button, thinner diameter grab bar	1"	4.69"	6.6"	5" W x 2.75"H x 15.2"L	100	1.7
	R1400229	Kid Grip fixed length 12.2", thinner diameter grab bar	1"	4.69"	12"	5" W x 2.75"H x 20.6"L	100	1.8
	R1400229S	Kid Grip 12.2" with suction indicator button, with thinner diameter grab bar	1"	4.69"	12"	5" W x 2.75"H x 20.6"L	100	1.9
	R1400555	Cardan joint adapter, pair, aluminum joint fits between pad and bar, swivel action moves in every direction, reduces bar capacity to 110 lb.		4.69"	Adds 5.25"	1" dia x 5.25"	110	.4
	R1400558	Hinged, flexible joint adapter, plastic joint with hinge fits between pad and grab bar, moves on one plane, useful in corners, Reduces capacity depending on length of bar.		4.69"	Adds 5.25"	1" dia x 5.25" L		.4
	R1400202	Shower head holder short, with suction pad. End of arm holds hand-held shower wand		4.69"		5" w x 2.75" H x 6"	8	.8
	R1400202S	Shower head holder short with suction pad and suction indicator button. Holds hand-held shower wand.		4.69"		5" w x 2.75" H X 6"L	8	.9
	R1499202	Shower head holder with swivel, suction pad with long arm swivels on pad, knuckle joint on end adjusts to hold hand-held shower wand	1.37"	4.69"		5" W x 7"H x 18"L	5	1.6
	R1499202S	Shower head holder swivel with indicator, suction pad with long arm swivels on pad, knuckle joint on end adjusts to hold hand-held shower wand	1.37"	4.69"		5" W x 7"H x 18"L	5	1.7
	R1400203	2 finger holder with suction pad, 2 articulated fingers can be molded to hold a variety of shapes		4.69"		5" w x 2.75" H x 11.15" L		.8
	R1400203S	2 finger holder with suction pad/swivel		4.69"		5" W x 2.75" H x 11.5"L		.8

	R1400204	Mobile clip holder		4.69"				1
	R1400204S	Mobile clip holder with indicator		4.69"				1.1
	R1499204	Mobile clip holder with swivel		4.69"				1.8
	R1499204S	Mobile clip holder/swivel with indicator		4.69"				1.9
	R1400205	Arm only with suction pad, without shower curtain		4.69"		5" x 11"W x 33"L		3.6
	R1400205S	Arm only with suction pad and suction indicator button, without curtain		4.69"		5" x 11"W x 33"L		3.7
	R14061208	Curtain only for R1400205				47" W x 31" H		
	R1400281U	Bathtub shortener, 2 suction pads	1.37"	4.69"		6" H x 16"W x 12"L	15	5.4
	R1400281SU	Bathtub shortener, 2 pads with suction indicator button	1.37"	4.69"		6" H x 16"W x 12"L	15	5.5
	R1400282U	Bathtub shortener, 3 suction pads	1.37"	4.69"		6" H x 16"W x 12"L	15	6
	R1400282SU	Bathtub shortener, 3 pads with suction indicator button	1.37"	4.69"		6" H x 16" W x 12" L	15	6.2
	R1400241U	QuattroPlus 7.9" - 12.8"	1.37"	4.69"	7.8" – 12.7"	14" – 18.9"	276	3.6
	R1400241SU	QuattroPlus 7.9" - 12.8" with indicator button	1.37"	4.69"	7.8" – 12.7"	14" – 18.9"	276	3.8
	R1400242U	QuattroPlus 12.6" - 17.5"	1.37"	4.69"	12.5" – 17.3"	18.7" – 24.4"	276	3.9
	R1400242SU	QuattroPlus 12.6" - 17.5" with indicator button	1.37"	4.69"	12.5" – 17.3"	18.7" – 24.4"	276	4.1
	R1400243U	QuattroPlus 16.6" - 21.5"	1.37"	4.69"	16.4" – 21.2"	22.6" – 27.5"	142	4.1
	R1400243SU	QuattroPlus 16.6" - 21.5" with indicator button	1.37"	4.69"	16.4" – 21.2"	22.6" – 27.5"	142	4.3
	R1400244U	QuattroPlus 25.4" - 30.3"	1.37"	4.69"	25.1" – 30"	31.4" – 36.2"	142	4.7
	R1400244SU	QuattroPlus 25.4" - 30.3" with indicator button	1.37"	4.69"	25.1" – 30"	31.4" – 36.2"	142	4.9

	R1400265U	QuattroPower support	1.37"	4.69"	22"	28.4"	176	6.5
	R1400265SU	QuattroPower support with indicator button	1.37"	4.69"	22"	28.4"	176	6.7
	R1400266U	QuattroPower support with extra grip	1.37"	4.69"	22"	31.2"	176	6.9
	R1400266SU	QuattroPower support extra grip with indicator button	1.37"	4.69"	22"	31.2"	176	7.1
	R1400261U	QuattroPower tub	1.37"	4.69"	33.9"	17"	176	6.2
	R1400261SU	QuattroPower tub with indicator button	1.37"	4.69"	33.9"	17"	176	6.4
	R1400262U	QuattroPower tub with extra grip	1.37"	4.69"	38.2"	17"	176	6.6
	R1400262SU	QuattroPower tub, extra grip w/ indicator button	1.37"	4.69"	38.2"	17"	176	6.8
	R1400206	Duo Cane holder, with suction pad and 2 rubber grips to hold canes		4.69"		5"W x 2.75" H x 10"L		1
	R1400206S	Duo Cane holder, with suction pad and suction indicator button. Has 2 rubber grips to hold a pair of canes		4.69"		5"W x 2.75" H x 10"L		1.1
	R1400209	Telescoping Mirror with suction pad, arm extends out 31.5". Mirror 6" dia. One side normal view, other side magnifies.		4.69"		2.75" H x 11"W x 31.5"L		1.1
	R1400209S	Telescoping Mirror, suction pad, with suction indicator button, arm extends out 31.5". Mirror 6" dia. One side normal view, other side magnifies.		4.69"		2.75" H x 11"W x 31.5"L		1.2
	R1409209	LED Micro light with clip, includes battery, clips onto mirror to direct a beam of light where needed						.1
	R1499209	Telescoping Mirror (details above) and LED light, with suction pad		4.69"		2.75" H x 11"W x 31.5"L		1.3
	R1499209S	Telescoping Mirror (details above), LED light, includes suction pad with suction indicator button		4.69"		2.75" H x 11"W x 31.5" L		1.3

	R14090222	Handi, stabilization handle, suctions onto table, vertical handle helps provide stabilization for hand tremors		3.75"		3.75" x 6"H		.5
	R1400295	Stabi horizontal stabilization handle, fixed handle on suction pad lays horizontal to flat surface, provides stabilization for hand tremors		4.69"		5" W x 2.75" H x 8.75"L		.9
	R1400295S	Stabi horizontal handle with indicator		4.69"		5" W x 2.75" H x 8.75"L		1
	R1400296U	Stabio Vario, adj. stabilization handle, handle is can be rotated 45 degree angles from horizontal to vertical from surface.		4.69"		5"W x 2.75" H x 10"L		
	R1400296SU	Stabio Vario, handle with suction pad and suction indicator button, handle can be rotated in 45 degree angles from horizontal to vertical from surface.		4.69"		5"W x 2.75" H x 10"L		
	R1461414	Stainless steel plates, set 2, provides the 5" x 5" flat surface needed for suction grab bars. Use areas with small tiles, grout lines, or porous surface. Adhesive backing				5" x 5"		