

# Clarke Health Care Products



## Bathlift Specifications/Comparison

	Aquatec	Aquatec RSB	Aquatec Wide	Aquatec RSB Wide	Aquatec XL	Major	Bali
Product number blue	A1471565	A1471575	A1471651	A1471582		A152112	
Product number white	A1471197	A1471202	A1471650	A1471579	A1471069		D10142
Length of seat, top of platform	19"	19"	19"	19"	19"	Custom	27"
Seat depth						Custom	
Width of seat	14.75"	14.75"	14.75"	14.75"	14.75"	Custom	11.75"
Width of top platform, opened	27.88"	27.88"	35.5"	35.5"	27.88"	Custom	
Height of back support	26"	26"	26"	26"	26"	NA	
Width of straight back support	14"		14"			NA	
Width of special back support, flaps open		27"		27"		NA	
Footprint, at bottom of tub	12" x 23"	12" x 23"	12" x 23"	12" x 23"	12" x 23"	Custom	
Maximum height	17"	17"	17"	17"	17"	Custom	
Height with adapters (1, 2 or 3")	18" - 20"	18" - 20"	18" - 20"	18" - 20"	NA	Custom	
Adjustable inner dimension							18.25" - 23.75"
Lifting capacity	300 lbs	300 lbs	300 lbs	300 lbs	400 lbs	up to 500 lbs	330 lbs
Weight of lift	32 lbs w/ battery	35 lbs w/battery	37 lbs w/battery	39 lbs w/battery	27 lbs w/battery	40 lbs - 60 lbs	4 lbs

## Accessories / Options

	Length	Width	Height	Weight	Thickness
Head Support for Major backs	19"	7.5"	11" - 16" above back	3 lbs	3"
Head Support, Beluga	8"	12"	8" - 10" above back	1.5 lbs	1"
Straight back for Major		12"	16"	3 lbs	.75"
Special back for Major		28.5"	16"	6 lbs	3"
Reclining straight back for Major		12.5"	17"	7.5 lbs	3"
Reclining Special back for Major		25.75"	17"	9 lbs	3"
Rotary seat with transfer board	21"	12.5"	21"	8 lbs	0.75"
Wedge Sand cushion	14"	9.5"	6"	18 lbs	
Pillow	14"	10"		1 lbs	2"
Seat belt	27"	2"			
Chest belt	31"	2"			
Battery	4"	4"	8"	3 lbs	