



The CORZO

Special people with special needs can travel in comfort and style. **Corzo** comes in 2 styles **Corzo Basic** in 3 sizes and **Corzo Cross Country (XCOUNTRY)** with lots of adjustability and 4 sizes. With a variety of sizes, **Corzo** can accommodate riders from child to young adult sizes. Optional accessories provide safe and comfortable positioning for these special travelers.

Features

- Folding steel frame
- Fashionable padded upholstery
- Split footrests
- Large tires for smooth ride
- Foot activated rear brake
- 5 point harness

Components

All **Corzo** models are built on a steel tubular frame that folds side to side, umbrella style. Hand grips are padded for caregiver comfort. The removable upholstery is padded and machine washable. The swing-away, flip-up, split foot plates can be adjusted in length. Calf strap attaches to foot plate with velcro. Deep traction wheels are 8" in front and 12" rear.

Corzo Basic includes the standard features above and has a Royal Blue frame and Black padded upholstery. **Corzo Basic** is available in 3 sizes 30, 34, 38. Size 38 has a detachable backrest extension to achieve the 33.5" height.

Corzo XCOUNTRY has all the above features and has growth capabilities with adjustments for back height, seat depth and push handles. **Corzo XCOUNTRY** has a Silver frame and Blue padded upholstery. **Corzo XCOUNTRY** is available in 4 sizes, 30, 34, 38 and 42. Tool for making height adjustments is stored in pouch on back fabric.

Warranty

Clarke Health Care Products and **Patron** provide a limited 2 year warranty against manufacturer's defects. Warranty does not include upholstery and tires.



SEAT ADJUSTMENT ON XCOUNTRY

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com





The CORZO Accessories



HEAD REST



SIDE LATERALS



FOOT STRAPS



ZIPPER HARNESS



H-HARNESS



FABRIC BASKET



HOOD

With this collection of accessories you can outfit the **Corzo** to your specific needs. All upholstered accessories are black hand-washable durable 100% Polyester. A retractable canopy hood shelters the head and face. Comfort seat cover provides extra padding. It will grow with the seat adjustments. Headrest can be adjusted for size, Foot straps attach to the foot plates and secure around the ankles with Velcro. Side lateral is a wide belt with built-in wedges to prevent side movement. H-style harness has quick release buckles at shoulder and hip. Zipper vest has a front zipper and shoulder and hip release buckles. Fabric basket attaches below the seat to carry personal items.

| SPECIFICATIONS | CORZO BASIC | | | CORZO XCOUNTRY | | | |
|-----------------------|-------------|--------|--------|----------------|--------|--------|-------------|
| | 30 | 34 | 38 | 30 | 34 | 38 | 42 |
| Backrest Height Min | 20" | 25.75" | 33.5" | 18" | 22" | 22" | 26.5" |
| Backrest Height Max | | | | 22" | 25.75" | 25.75" | 30.5" |
| Seat depth Min | 10" | 15" | 15" | 10" | 11.75" | 12.5" | 13.75" |
| Seat depth Max | | | | 13.25" | 15" | 15.5" | 15.5" |
| Seat width | 11.5" | 13.25" | 15" | 11.5" | 13.25" | 15" | 16.5" |
| Seat to footrest Min | 6.5" | 10.5" | 10.5" | 6.5" | 10.5" | 10.5" | 10.5" |
| Seat to footrest Max | | | | 14.25" | 18.25" | 18.25" | 18.25" |
| Backrest angle | 115° | 115° | 115° | 115° | 115° | 115° | 115° |
| Seat-unit angle | 25° | 25° | 25° | 25° | 25° | 25° | 25° |
| Total length Min | 43.25" | 44.5" | 44.5" | 44" | 44.5" | 44.5" | 44.5" |
| Total length Max | | | | 40.25" | 40.75" | 40.75" | 40.75" |
| total width | 21.5" | 23" | 24.75" | 23" | 24.5" | 26" | 38.5" |
| Total height w/handle | 36" | 40" | 40" | 37.5" | 38" | 38" | 30"- 40.75" |
| Weight (lbs) | 26.5 | 27.5 | 31 | 25 | 26.5 | 30 | 31 |
| Capacity (lbs) | 100 | 120 | 165 | 100 | 120 | 165 | 165 |
| Folded length | 45" | 48" | 48" | 45" | 46" | 46" | 47" |
| Folded width | 14" | 14" | 14" | 15" | 15" | 15" | 15" |
| Folded height | 12" | 14" | 14" | 14" | 14.5" | 13" | 13" |

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

