

Arcorail™ rotating 360° bedside handrail professional grade



The rotating, professional-grade bedside handrail gives you maximum bedside mobility for home, clinical care and hospital settings. Arcorail™ fits securely to the metal bed frame and helps you reposition, sit-up, balance, stand or transfer to a chair or wheelchair.

- Rotates 360° & locks in place every 90°
- Professional grade; made of steel & ideal for frequent use
- Special coating does not promote bacteria growth
- 3 height adjustable settings
- 300 lbs max. user weight capacity


ADJUSTABLE FOR USE WITH HOME STYLE OR CLINICAL CARE & HOSPITAL BEDS

HOME STYLE BEDS

Fits angle-iron bed frame used with box spring and mattress (twin to king, standard to pillow-top)

CLINICAL CARE & HOSPITAL STYLE BEDS

Fits Clinical Care style beds with inverted angle-iron bed frame where the head and foot are adjustable and Hospital style (long-term care) beds with round or square tubing frame, such as beds made by Hill-Rom.



U.S. Patent 5,384,927 / Canadian Patent 2,096,681
U.S. Trademark 1,860,819 / Canadian Trademark 518,317

www.hartmobility.com

Hartmobility™
Premium HOME HEALTHCARE

New Compact Packaging



Box Dimensions = (W) 14 1/2" x (H) 22 1/2" x (D) 4 3/4"



CLARKE
HEALTH CARE
PRODUCTS INC.

Toll Free: (888) 347- 4537
info@clarkehealthcare.com

Hartmobility™
Premium HOME HEALTHCARE
www.hartmobility.com