



# CARE BAG

*What can you say about bedpans, they stink and they are disgusting. No one wants the job of emptying the pan or pail. But you never have to spill a commode pail again. Simply line the collection pail with a **Classic Care Bag** and cleanup is easy. The **Care Bag** system eliminates spilling, odor and the spread of bacteria caused by handling a full bedpan or pail.*

### **Don't have a bedpan?...**

*Use a VOM bag for upset stomachs.*

*URI bags are used as a personal male urinal.*

### **Features**

**Care Bags** come in sizes to fit commode pans and pails, or for personal use. Each bag comes with a GelMax absorbent pad that quickly absorbs up to 16 oz. of human liquids and odors. All liquid collected is converted to a gel. Cleanup is as easy as tying the bag closed and disposing in the trash.

- 3 Styles
- Eliminates or reduces odors
- Liquid is converted to gel, no spilling
- 30 second bedpan cleanup



*Completely line/cover a bedpan or collection pail.*



*Liquid is absorbed and converted to gel.  
(Color added to enhance photograph.)*

**CLARKE HEALTH CARE PRODUCTS, INC.**

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

CARE BAG



# CARE BAG



CARE BAG CLASSIC  
C7711141

CARE BAG CLASSIC S  
C7831738

CARE BAG URI  
C7733993

CARE BAG VOM  
C7733987

Hospital staff or caregivers can avoid the spread of infection caused by spillage or contact with the body fluids. With the recent attention to hospital acquired infections, the value of **Care Bags** can be seen in decreased staff downtime due to HAI and improved patient safety. The use of **Care Bags** does not require training and can easily fit into an existing hygiene procedure.

## Components

**Care Bag** systems are polyethylene bags with pull tie closures. Each bag comes with a GelMax absorbent pad. Double-soldered seams along with the pad means no risk of leaking. Classic and Classic S bags are sized to fit most commode pans, bedpans or pails. URI bags are sized to collect up to 16oz of urine discreetly. VOM bags are sized to be placed over the mouth to collect vomit.

## Specifications

PRODUCT NUMBER	DESCRIPTION	SIZE	WEIGHT
C7711141	Classic Care Bag	24" x 17"	1.4 oz ea / 1.3lb pk
C7831738	Classic S Care Bag	19.25" x 14"	1.2 oz ea / 1.5 lb pk
C7733987	VOM Care Bag	9.5" x 9.5"	.7 oz ea, / 1 lb pk
C7733993	URI Care Bag	5.25" x 12"	.8 oz ea / 1 lb pk

## CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071  
Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

CARE BAG