

# Portable aids

**MOBELI®** Cane holder • **MOBELI®** Clip holder

The times are past when your cane or crutches were lying on the floor or leaning in some corner at home. The **MOBELI®** Duo cane holder with particularly strong adhesive strength can be fixed on a flat, non-porous wall surface and also on or under a table very easily. The ROTH-**MOBELI®** Clip holder with 2 Finger holder – Cane holder or Shower holder will be clipped on the **MOBELI®** Grips page 2 + 3.

<b>MOBELI®</b>	Order no.
8 Duo-cane holder	1400 52-90
9 Clip-holder for cane*	1400 62-90
10 Clip-2 finger holder for hairdryer e.g.*	1400 62-93
Clip-shower holder (see picture above)*	1400 62-02

\* EP angem.



7

**MOBELI®** Telescope-mirror  
**MOBELI®** Shower holder  
**MOBELI®** 2 finger holder



Just as practical aids are the **MOBELI®** 2 finger holder and the **MOBELI®** shower holder – The **MOBELI®** Telescope-mirror will be a good help for Stoma carrier and with the LED light it will be a help to catheterize. Each product will be able to rotate and be fixed in the optimal position.

<b>MOBELI®</b>	Order no.
11 MOBELI Telescope Mirror	1400 52-09
12 LED-Micro Light with clip for Telescope mirror	1409 0209
13 Two Finger holder e.g. for hairdryer	1400 52-93
14 Shower head holder long*	1400 52-03
15 Shower head holder short*	1400 52-02

\* EP angem.



8

# Stabilization handles

**MOBELI®** Stabi • **MOBELI®** Handi • **MOBELI®** Kid-Therapy-Grip "Daisy"

Wherever you are, the practical **MOBELI®** Handi (picture 18) offers fast help for immobilizing your arm/hand, e.g.: with Parkinson, ataxia, spastics, etc. **MOBELI®** Handi can be simply attached to any smooth and non porous area, and just as easily removed again. A further benefit: With 80 mm diameter and its low weight, **MOBELI®** Handi is also ideal for traveling. The detachable stabilization grip **MOBELI®** Stabi is particularly versatile in its application. In contrast to the simple horizontal grip part (picture 17), the Vario (picture 16) grip part can alternatively be attached in five different positions and thus adapted to the respective requirements. On a smooth area, the 120 mm suction plate develops twice as much adhesive strength as its smaller brother **MOBELI®** Handi.



9

The **MOBELI®** Kid-Therapy-Grip "Daisy" Vario (picture 19) also has five different positions and with the "Daisy" horn will promote a child's fine motor skills. The children will squeeze the blue hat or the orange beak of "Daisy", it will help them when they are in motion. The "Daisy" horizontal (picture 20) has a grip length of 160 mm, 6.3" and can be gripped by 2 small hands.

<b>MOBELI®</b>	Grip-Ø in inch	Adhesive strength *	Order no.
16 Stabi Vario	1.2	200 N (~20 kg, 44 lbs)	1400 52-98
17 Stabi Horizontal	1.2	200 N (~20 kg, 44 lbs)	1400 52-95
18 Handi	1.2	100 N (~10 kg, 22 lbs)	1409 0222
19 Kid-Therapy-Grip Daisy Vario**	1	150 N (~15 kg, 33 lbs)	1400 52-88
20 Kid-Therapy-Grip Daisy Horizontal**	1	100 N (~10 kg, 22 lbs)	1400 52-89

\*\* EP angem.



10

# Mobile hand grips

For at home and travelling.



**MOBELI®**  
**MOBILITY**  
**BEYOND**  
**LIMITS**

**ROTH®**  
Sanitärprogramme  
REHA products

# Mobile hand grips

**MOBELI®** Single and dual hand grips • **MOBELI®** Telescope • **MOBELI®** Kid

Control your everyday life!

How fast does something happen – at sports, at home, at work – that afterwards impedes us in our life rhythm for weeks, for months, or a lifetime.

The **MOBELI®** system by ROTH is a practical aid that has been proven in the rehabilitation and in everyday life of handicapped or older persons in the best way for many years.

The mobile handles are quickly attached – **without drilling** – and just as easily released again. The two suction plates with 4.6" (120 mm) diameter, are made of a special rubber mixture and are simply pressed onto any smooth, porous surface. The vacuum toggle lever is pushed down – that's it! The handles will hold with high suction strength and offer support wherever you need it, both at home and when traveling. Always test the holding grip on each surface, each time you apply the suction pad.

	<b>MOBELI®</b>	grab length x Ø inch/mm	Adhesive strength *	Order no.
1	Single grip	4"/100 mm x 1,4"/35mm	700 N (~70 kg, 154 lbs)	1403 999
1	Dual grip	13.25"/340 mm x 1,4"/35mm	600 N (~60 kg, 134 lbs)	1400 52-6
1	Dual grip	25.2"/640 mm x 1,4"/35mm	600 N (~60 kg, 134 lbs)	1400 52-7



2


For you, our claim of quality means security that you can trust: Our handles have been awarded the TÜV seal.

You can disassemble the mobile dual grips to carry with you when traveling: The suction plates are easily removed from the grip bar without tools.

The mobile hand grips are available in selected lengths, design variations, and with practical accessories.



	<b>MOBELI®</b>	grab length inch/mm	Adhesive strength *	Order no.
2	Dual grip telescope	8.25–13.2"/210–335	600 N (~60 kg, 134 lbs)	1400 52-1
2	Dual grip telescope	13–17.9"/330–455	600 N (~60 kg, 134 lbs)	1400 52-2
2	Dual grip telescope	16.9–21.8"/430–555	600 N (~60 kg, 134 lbs)	1400 52-3
	Kid grip	8–8.9"/205–225 x Ø 1"/26 mm	600 N (~60 kg, 134 lbs)	1400 52-8



3

# Portable Bath assistance

**MOBELI®** QUATTRO+PLUS • **MOBELI®** Bath tub shorteners • **MOBELI®** Universal joint

The ROTH-**MOBELI®** grip QUATTRO+PLUS\* with 4 suction plates, 120 mm (4.6"), has the integrated patented ROTH®-Universal Joint and 100–120 kg (220–260 lbs) holding power. The grip can be used on the wall or in a corner. It is possible to use the grip as a training grip between a pore-tight metal door frame. The grip length of 645–770 mm (25.4–30.3") will fit a door frame of 790–905 mm (31.1–35.6").

	<b>MOBELI®</b>	grab length inch/mm	Adhesive strength *	Order no.
	QUATTRO+PLUS**	7.9–12.8"/200–325	1200 N (~120 kg, 260 lbs)	1400 52-1Q
	QUATTRO+PLUS**	12.6–17.5"/320–445	1200 N (~120 kg, 260 lbs)	1400 52-2Q
4	QUATTRO+PLUS**	16.5–21.5"/420–545	1000 N (~100 kg, 220 lbs)	1400 52-3Q
3	QUATTRO+PLUS**	25.4–30.3"/645–770	1000 N (~100 kg, 220 lbs)	1400 52-4Q

\*\* US Patent pending



4



The application range of the **MOBELI®** telescope handles can be substantially extended by universal joint adapters which are installed between suction plates and telescope bar if necessary – the adhesive strength will be reduced to 40 kg (88 lbs). With the installed universal joint the **MOBELI®** telescope handles can be used at the corner and very universally mounted. The infinitely variable bath tube shortener which is also be usable as a tray (see picture above) can be used for most bathtubs with 2 or 3 (for higher power) suction devices.

	<b>MOBELI®</b>	Tub width/ inch/mm	Suction devices	Order no.
5	Bath tub shortener/tray	20–27.6"/510–700	3	1400 52-07
	Bath tub shortener/tray	20"/510	2	1400 52-06
6	Pair of universal joints			1400 555



5

You can find additional ROTH® products by calling Clarke Health Care. These are special order items:



Clarke Health Care  
1003 International Drive  
Phone 888 347 4537  
Fax 724 695 2922  
info@clarkehealthcare.com  
www.clarkehealthcare.com

Quick grip and hold ... Quick fitting for bathtub  
Now you can install a grab rail directly onto your bathtub to help you get into the bath. Also you can use it for the corner. With the telescope version you can use it for various tiles. The ROTH®-Patent Combigrip Mobile Bath Grab Rail can be converted to a mobile holding grip without needing any tools – see the information page 2 + 3.

The bath grab rail with its 4 suction devices can support 134 lbs/ 60 kg if used appropriately.

	<b>MOBELI®</b>	grab length inch/mm	Order no.
	Combigrip Bath Grab Rail*	11.8"/300	1400 52-6W
7	Combigrip Bath Grab Rail*	11.4–16.4"/290–415	1400 52-2W
	Combigrip Bath Grab Rail*	15.4–20.3"/390–515	1400 52-3W

\* US Patent pending



10/06 Änderungen vorbehalten.

D-72213 Altensteig  
Phone 00 49 74 53 / 93 81-0  
Fax 00 49 74 53 / 93 81-23  
info@roth-reha.de  
www.roth-reha.de

**ROTH®**  
Sanitärprogramme  
REHA products

6