


Clarke Health Care Products

Shower Chair Specifications/Comparison

								
	Ocean	OceanXL	Ocean SP	Ocean VIP	OceanEVIP	Dukki2	Dukki3	Dukki 4
Product number	A710900	A713900	A712900	A1470713	A1471208	H21080	H31080	H41080
Back height	14"	14"	14"	21"		16.5"	18"	19"
Seat height	19"-23"	19"-23"	19"-23"	19.5"/21.5"/25.75"	19.5"-41"	19-24"	18-23"	18-23"
Seat width	19.5"	19.5"	19.5"	19.5"	18.75"	11"	12"	13"
Seat depth	18"	18"	18"	17.5"	17.5"	11"	13"	15"
Overall width	22.5"	24.5"	27"	22.5"	28.5"	22"	23"	23.5"
Overall length	37"	37"	42.5"	41.25"	39.5"	23"	26"	28"
Dist seat to footrest	12" - 22"	12" - 22"	12" - 22"	14"-19"		9-12"	10-15"	11.5-17"
Hght under seat w/o pan	17.25" - 22.5"	17.25" - 22.5"	17.25" - 22.5"	16.5" - 20.25"		17.5" - 22.5"	16.5" - 21.5"	16.5" - 21.5"
Width between armrests	18"	21"	18"	17.5"	18"	NA	NA	NA
Dist between rear wheels	15"	15"	15"	18"	17.5"	17"	18"	18.5"
Tilt angle forward	NA	NA	NA	0°		-5°	-5°	-5°
Tilt angle backward	NA	NA	NA	35°		-60°	-60°	-60°
Max weight capacity	290 lbs	330 lbs	290 lbs	330 lbs	330 lbs.	80 lbs	80 lbs	80 lbs
Weight	33.5 lbs	33.5 lbs	47 lbs	48 lbs	83.75 lbs.	22 lbs	26.5 lbs	29 lbs

Product name	Product number	LxWxH	Weight	Fits chair	Notes
Reclining back support	A1470085			OceanVIP	will not fit EVIP
XL back support	A16515	19" x 21" x 18"	7.5 lbs	Ocean, OceanSP	adds 3" arm width
Safety Harness combination	D10068	1" wide	1 lbs	all Oceans	set of 3 belts
Safety Support bar	D10167	22" x 1" dia	2 lbs	does not fit XL	
Calf strap	D10175	9.2" x 4" x 1.25"	.5 lbs	all Oceans	Velcro closure
Hygiene Seat Overlay	A1470075	18" x 16" x 3.25"	5 lbs	OceanVIP, OceanEVIP	
Hygiene Seat Overlay	A16341	17.5" x 17.5" x 3"	4 lbs	Ocean, SP, XL	
Toilet Seat Overlay	A16342	17.5" x 17.5" x 3"	4 lbs	Ocean, SP, XL	cutout 9" x 6.25"
Action Gel	D100031	18" x 18"	8 lbs	all Oceans	cutout 10" x 7"
Splash guard	A16511	4.5" x 4.5" x 10"	1 lbs	only hygiene seat	
Solid Seat Insert	A16343	15" x 8" x 2.5"	2 lbs	only hygiene seat	
24" wheels, pair	A16379		12 lbs	Ocean or OceanXL	converts to SP

Carry bag	A16521	16" x 20" x 25"		does not fit VIP, EVIP	
Padded chest belt	A1470081	4" x 30"/34.5"/38"	1 lbs	all Oceans	
Padded lap belt	A1470082	4" x 24"/29"/34"	1 lbs	all Oceans	
Head support		9" x 4"		included with VIP/EVIP	
Footpad		6" x 5.25"		included all chairs	