



# The LEGACY Valor

Walk with pride and ease, **Valor** is based on the popular **Dolomite Legacy** model but with extended handles, wider seat and Low-Vision paint scheme. The handles will adjust from 30" high to 39" high. A wider frame provides a more comfortable seating area. The sporty red color frame is complimented by the black seat and black spoke wheels.

## Components

The strongest frame on the market, **Dolomite** welds the steel frame for greater durability, no hardware at the joints to loosen. A corrosion-resistant finish is applied and then a powder-coat color is baked on. **Valor's** frame is easily folded for storage or transportation. The aluminum handle tubes have an 8" infinite adjustment range with lines marked along the tubes so you can duplicate the position.

The large seat flips up for easier access. Push the seat down to hear an audible snap assuring the frame is now locked in the open position; a red safety latch must then be released to close the walker. Large 7 1/2" wheels ride easily over varied terrain both indoors and outdoors. A flexible, adjustable back strap is included. The wire basket can hold personal belongings, be removed and carried with its own bail handle.

**Valor** uses the patented **Dolomite** Automatic Braking System with adjustable tension for varied grip strength. Brake handles have wide flanges and large loops for comfortable hand positioning. Push down the parking brake when parked and seated. The brake pad is located within the rear forks for protection.

## Features

**Valor** supports 300 pounds. Other standard features include:

- Flip-up seat for greater reach or gait training
- Seat handle release to fold walker
- Safety latch prevents accidental closure
- Foam grips for comfortable and secure hand position
- Removable basket with carry handle



Handle adjustment



## CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)



# The LEGACY Valor

## Accessories



OXYGEN TANK HOLDER



BAG



CANE HOLDER



ONE-HAND  
BRAKE



ANATOMIC  
HAND GRIPS



SLOW DOWN BRAKE



FOOD AND BEVERAGE  
TRAY



BASKET

## Warranty

Dolomite offers a limited warranty against manufacturer's defects; the frame is covered with a limited lifetime warranty. The DABS brake handle is covered for 5 years, other parts are covered for 1 year. This warranty excludes brake cable, foam grips and wheels.

## Specifications

SPECIFICATIONS	LEGACY VALOR
Handle Height	30" - 39"
Width Between Handles	18"
Seat Dimensions (L x W)	17" x 8.5"
Seat Height, from floor to seat	24"
Distance Between Rear Wheels	19.5"
Overall Width x Length	24.5" x 27"
Folded Dimension (L x W x H)	32" x 24.5" x 11"
Walker Weight	16.5 lbs
Weight Capacity	300 lbs
Product Number	D1522398

## CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

