

AQUATEC® *The* BALI

Comfortable and stable bathboard.

Simple and inexpensive alternative to bath chairs, **Bali** bath board sits across the top of your tub. Inside adjustable grips fasten the seat to the edge of the tub. Non-slip, textured seat gives you a secure surface for bathing. **Bali** has rounded edges for snag-free transfers. Hand grip offers additional safety. Recess on end of seat holds soap or bottles.

Bali is made from molded plastic components, white color. Simply place on top of tub and tighten the grips under the seat. Clean surface with ordinary bathroom cleaners.



NON-SLIP SURFACE

Warranty

Aquatec warrants this product to be free of manufacturing defects for 2 years.

SPECIFICATIONS

Part Number	A10142
Total length	27"
Total width	11 3/4"
Adjustable inner dimension, (inside top opening of tub)	18 1/4" - 23 3/4"
Width of hand grip handle	8 1/2"
Weight of Bathboard	4 lbs.
Weight capacity	330 lbs.

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

The BALI