



# The BANJO

Take a walk with a **Banjo** and enjoy its modern design. A substantial frame, with individual adjustments, provides the customization you need when walking. With the stand-up storage feature, **Banjo** will definitely stand out in the crowd.

Your **Banjo** rollator has great cargo space. It has a plastic shelf to hold up to 15 pounds of bulky items. A wire basket hangs from under the seat, easily detached and carried by the handle. A plastic tray secures onto the seat, the deep center well will contain spills.

Anatomic grips are comfortable and oversize loops on the hand brake handles provide large surface to activate the brakes. The extended adjustment on the handle height gives you a greater range of personal comfort. When you tire, sit and rest on the **Banjo** seat. Be confident that the wide tires will take you smoothly over varying floor surfaces or outdoor terrain.

## Features

- Steel frame
- Large indoor-outdoor wheels
- Detachable shelf, tray and basket
- Handles adjust over 7"
- 287 pound capacity

## Components

**Banjo** is crafted with a steel frame and chrome handle bars. There is a 7-inch range of handle height adjustment, secured with a locking lever. Formed anatomical grips and rounded brake handles provide good grip when braking. A redesigned safety catch makes it easy to release and fold the rollator. Large 7 1/2-inch tires for indoor or outdoor travel. A plastic shelf has 15-pound capacity, wire basket with hinged bail handle, impact-resistant plastic tray attaches securely on the seat.



**CLARKE HEALTH CARE PRODUCTS, INC.**

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)





# The BANJO Accessory

## Accessory

An optional cane holder clips onto the frame. Two-piece design, the bottom cup will cradle the foot of your cane.

## Warranty

**Banjo** is covered by a two year warranty against manufacturer's defects. This warranty excludes wear and tear parts such as wheels, hand grips and brake cables.



*Banjo shown with optional cane holder*

## Specifications:

SPECIFICATION	BANJO, P452E-3
Overall width	23.5"
Overall length	27"
Width between rear wheels	17"
Height	27.5" x 38.25"
Seat	15"L x 6.5"W x 1"H
Seat height from floor	23.75"
Wheel	7.5"D x 2"W
Weight	19.75 lbs
Capacity	287 lbs
Shelf capacity	15 lbs
Shelf size	16.75"L x 8"W x 1"D

## CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • [www.clarkehealthcare.com](http://www.clarkehealthcare.com)

