

AQUATEC® *The* OCEAN VIP

The Ocean VIP with New Recline Back Option

The next wave in comfort and style, **Aquatec's Ocean VIP** tilting shower chair is height adjustable, for easier transfers or toilet clearance. **Ocean VIP** provides correct posture for your hygiene needs with tilt-in-space or recline positioning choices. Recline positioning now available with the optional back support.

Features

- Stainless steel frame
- Seat height easily adjusted from 19.5" to 25.5" from floor
- Adjustable headrest
- Adjustable, machine washable back upholstery
- 5" casters with locks
- Height adjustable, flip-up armrests

Components

Ocean VIP is built on a stainless-steel frame with corrosion-free components. The seat is waterfall style with a hygiene cut-out for comfort. Soft overlays give you seating choices for toileting or showering. Two stainless-steel gas springs will lift 330 pounds. Adjustable seat height is set with easy-to-move clips. Footrests can be lengthened easily with locking pins and include tilting footplate and heel loops.

The back upholstery has 5 buckles to adjust tension for seating comfort. Optional reclining back support provides even greater seating and positioning possibilities. The angle between the reclining backrest and seat is infinitely adjustable from 10° to 50° for a total angle of 85° for the complete chair. The operating lever locks to avoid unintentional angle movement. A pair of anti-tippers is included with this back option.

The **Ocean VIP** rolls smoothly on 5" casters. Three wheels have a foot activated brake and one rear wheel has direction lock for straight line travel. A handle on the back push bar sets the tilt from 0°-35°. The gas springs are stainless steel for long service. Armrests flip-up for side transfers and can be positioned at two heights. Optional arm locks are available.

Warranty

Ocean VIP shower commode chair has a limited 2 year warranty against manufacturer's defects. The optional reclining back support has a limited 2 year warranty against manufacturer's defects.



*Seat Height
Adjustment*



*Armrest
Adjustment*



CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

AQUATEC® *The* OCEAN VIP



TOILET OVERLAY



SOLID INSERT



SPLASH GUARD



GEL SEAT



PADDED SEAT AND CHEST BELT



SAFETY SUPPORT BAR



SAFETY HARNESS COMBINATION



ARM LOCK



ARM WIDTH EXTENSION



SIDE LATERAL AND ARM TROUGH



CALF OR AMP SUPPORT



CALF STRAP



CARRY BAG

SPECIFICATIONS

	OCEAN VIP
Back height from seat	21"
Seat height from floor	19.5"/21.5"/25.75"
Seat width	19.5"
Seat depth	17.5"
Overall width of chair	22.5"
Overall length 0° tilt	41.25"
Overall length, full tilt, footrest at max, headrest at max	60.2"
Width between arms	17.5"
Height floor to lower seat edge w/o pan holder	17.25" - 22.25"
Tilt angle	0°-35°
Distance between rear wheels	18"
Max weight capacity	330 lbs.
Weight	48 lbs.
Product Number	A1470713

RECLINING BACK SPECIFICATIONS

Width between arms at end of handgrip	18.5"
Reclining back angle	0°-50°
Overall width of back	22.5"
Overall length full tilt, no recline push bar to foot plate (no head support)	54"
Overall length, full tilt and recline, head support to foot plate	68"
Height from seat to top of upholstery	20"
Length chair with recline back, no tilt, no recline, from anti-tipper to footrest	42"
Product Number	A1470085

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com