



The PARAPLIS Umbrella Holder

Rainy days are no longer a problem when you're out with a stroller or wheelchair. **Paraplis** attaches to your stroller or wheelchair to hold an umbrella beside you as you walk. When in use, it swings up to hold an open umbrella. When not in use, it swings down to carry the closed umbrella at the side of the stroller.

Features

- Attaches to any shape bar or tube
- Adjustable clamp holds open umbrella upright
- Easy push button release swings holder up or down
- Rubber pads cushion the tube or umbrella handle
- No loose parts to get lost



Components

Paraplis consists of a 7 1/2" reinforced plastic bar with two adjustable clamps. Each clamp has a push button release for angle adjustments. Rubber inserts inside each clamp prevent damage to the wheelchair or umbrella. Easy thumbscrew hinges into place to secure the clamp around the wheelchair tube or the umbrella stem. Clamps will fit tubes from 3/4" – 1 3/8" diameter. There are no loose parts to get lost. Bottom of umbrella handle will extend 4" above bottom clamp.

Warranty

Paraplis is covered by a one year warranty against manufacturer defects. Caution, **Paraplis** is not intended to hold an open umbrella in high wind situations.

Specifications

SPECIFICATIONS	PARAPLIS
Total Length	7 1/2"
Total Width	5"
Weight	9 oz.
Fits tube/umbrella diameter	3/4 – 1 3/8"
Product Number	EPARAPLIS081



CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com