



The SOPRANO

Amid a chorus of Ohh's and Ahh's, the **Soprano** has joined the line of beautifully engineered **Dolomite** walkers. **Soprano**, with its larger wheelbase, provides a smooth ride over life's bumps. Completely redesigned with new frame, seat, handles and cable-less braking system.

Frame

Soprano is based on a welded steel frame. The frame is designed to eliminate the rear forks, giving you more walking space between back wheels. A larger seat flips up to give you more walking room. **Soprano** folds in two positions, first to stand upright by itself and then a second position to latch securely for travel. The front wire platform can carry up to 26 lbs of baggage.

Braking

The brake handles and wedge shaped hand pad offer a secure and comfortable grip. Cable-less braking system uses stainless-steel rods encased in the frame to actuate the brake. The brake shoe is designed to catch dirt and stones before they get stuck in the brake pads, adding longevity to the wheels and brakes. Infinitely height-adjustable handles also rotate for better wrist positioning.

Standard Features

- Welded steel frame
- 10" wheels
- Wire shelf carries your bulky items
- Handles rotate and are height adjustable
- Stands when folded, latches for storage
- Pearl Blue-Graphite frame color

Warranty

Soprano is guaranteed for five years.



STANDS IN
FOLDED POSITION

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com





The SOPRANO

Accessories



BASKET



CANE HOLDER



WHEEL COVER



CURB CLIMBER



TRAY



O2 TANK HOLDER



SLOW DOWN BRAKES

The following accessories have been developed to enhance your **Soprano**. Hold smaller items with the accessory basket. A cane holder attaches to the frame to hold a cane or crutch. The wheel cover keeps dirt or mud from splattering on your clothing. Use the curb climber to assist when tilting the walker to climb over curbs. Food and beverage tray sits on the seat, will also hold small items. You can safely carry a B, C or D size tank with the Oxygen Tank Holder. Slow-down brakes add tension to the wheels so the walker won't roll away, useful for gait training or hilly terrain. The tension is self-adjustable and can be secured in position.

SPECIFICATIONS	SOPRANO
Inside width between rear wheels	20"
Min. handle height	30.5"
Max handle height	37"
Seat height	23.5"
Width between handles at seat	17"
Overall length of walker, front to back	33.5"
Overall width	25.5"
Folded dimensions	11.5"L x 25.5" W x 34.5" H
Walker capacity	355 lbs.
Walker weight	21.4 lbs.
Product number	D12030

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

