



The TOM

Pediatric Positioning Transportation Systems

Cruise town and country in style. **TOM** strollers and accessories offer a variety of seating possibilities for the young rider. Available in 4 sizes, *Mini, Standard, Maxi and Super Maxi*, **TOM** strollers will fit children and young adults up to 165 pounds. Tilt or recline in comfort on fashionable padded upholstery. The seating system tilts or reclines. **TOM** strollers not only provide a comfortable ride but also convenience for the caregiver, with an ergonomic adjustable handle bar.

Components

TOM strollers have a removable seat system on a steel chassis in a Silver powder-coat finish. Green/Charcoal Tweed padded upholstery is machine washable. An upholstered grab bar is included. The padded push handle bar is adjustable on 3 models. The spring-suspension chassis has 12" rear wheels with deep traction for outdoor travel. Rear foot-activated brake locks the rear wheels. The 8" front casters have a locking swivel.

Features

The removable seat system can be reversed for more interaction with the child. Frame folds for storage or stowing in your vehicle and can be folded with the seat and accessories in place. The seating system has a reclining or tilt-in-space movement. **TOM** has passed EN ISO 7176/19 crash test standards and includes transportation tie-downs.

TOM frame and seat have multiple adjustment points.

- Height-adjustable back support
- Back goes to full recline
- Seat tilts +/- 27°
- Reversible seat direction
- Seat depth adjusts 3"
- Seat width adjusts up to 3"
- Footrest length adjusts up to 6"



TOM folded

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com



The TOM

Accessories

Optional accessories will make your rider very comfortable and stable. They include head supports, side supports and pads, foot straps, a variety of belts and harnesses, hoods and covers, basket and bags. Upholstered accessories color coordinate with the seat upholstery. For complete details see separate sheet.



Warranty

Clarke Health Care Products and Patron warrant this product to be free of manufacturing defects for 2 years, excluding wheels and upholstery.

SPECIFICATIONS	TOM MINI	TOM STD.	TOM MAXI	SUPER MAXI
Backrest height	17" - 23.75"	20" - 25"	22" - 26"	22" - 26"
Seat depth	9.5" - 12"	9.5" - 12"	12.5" - 15.5"	12.5" - 16.75"
Seat width	9.5" - 11"	9.5" - 11" - 12.5"	11" - 12.5" - 14.5"	13" - 15" - 16.25"
Seat to footrest	5" - 12"	7.75" - 13.75"	9.25" - 15.5"	8.5" - 15"
Backrest angle	0 - 180°	0 - 180°	0 - 180°	0 - 160°
Seat angle	(+/- 27°)	(+/- 27°)	(+/- 27°)	(+/- 25)
Total length	33.5"	33.5"	36"	43"
Total Width	24"	24"	25.5"	27.5"
Total Height/ with handle	39" - 45"	39" - 45"	40" - 46"	42"
Weight (lbs)	44	46.3	52	53
Weight capacity (lbs)	77	77	110	165
Folded stroller length	36.5"	36.5"	37.5"	37.5"
Folded stroller width	24"	24"	25.25"	27.25"
Folded stroller height	21.5"	21.5"	24.5"	25.5"

CLARKE HEALTH CARE PRODUCTS, INC.

1003 International Drive • Oakdale, PA 15071

Phone (888) 347-4537 • Fax (724) 695-2922 • www.clarkehealthcare.com

