

Safety in the spotlight

MOBELI® - always a thought through concept.

The innovative original.

MOBELI®

Position. Stick. That's it.
Safety in the spotlight.



Clarke Health Care Products, Inc
1003 International Drive
Oakdale PA 15071
Toll Free: 1-888-347-4537
www.clarkehealthcare.com

ROTH®

D-72213 Altensteig
Telephon +49 (0) 7453 9381-0
Fax +49 (0) 7453 9381-23

www.mobeli.de

You can reach us at info@mobeli.de

Mobile helpers

MOBELI® - small helpers, big impact.

Our mobile helpers support you in a variety of ways in daily living. Just position our helpers in the desired place to mount a hairdryer, shaver or hair-brush – that's it.

The helpful shower positioner is available with a flexible swivel arm as desired. The mobile telescopic mirrors – with optional LED light – are fast and easily positioned. They make it possible to look into the mirror while seated or are very helpful during catheterizing.

If your canes are on the floor too often, the mobile Duo cane holder will make life easier for you. It can be mounted on the wall or under the tabletop. It's always fast and safe and you never need tools.



MOBELI®	Order no.
Mobile clip holder with swivel	14992 04 S
Duo cane holder	14002 06 S
Shower head positioner with swivel	14992 02 S
Shower head positioner	14002 02 S
Mobile telescope-mirror (31.5")	14002 09 S
LED micro light with clip (incl. 50 h battery)	14092 09



The name says it all: visual screen and splash-guard; very popular and ideal in nursing care. The hinged bracket can be swivelled upwards or can be panned to the side when vertically mounted. The integrated brake holds the Mobeli arm in the desired position automatically. The arm was designed in a way that weight applied puts less pressure on the sustaining wall due to its long three-part elastic design.



MOBELI®	Order no.
Arm for protective curtain (without curtain)	14002 05 S
Antibacterial curtain (47.2" x 31.5")	140612 08

Support systems

Safer grip. Exactly where you need it.

Give yourself instantly a safe grip without drilling. Position the QuattroPower support system, press down the rocker, that's it. Now you have a reliable aid on flat and non-porous surfaces.

The QuattroPower Support system with 4 vacuum pads and a very special construction for reducing the leverage can carry a nominal load of up to 176 lbs and can be used even on tiled hollow space walls.



MOBELI®	Order no.
QuattroPower Support	14002 65 S
QuattroPower Support with grip	14002 66 S
QuattroPower Tub	14002 61 S
QuattroPower Tap with grip extension	14002 62 S



Optimized for the bathtub: The QuattroPower tub is perfect for getting in and out of the tub. It can sustain pull or bearing loads of up to 110 lbs nominal load.

Did you know that Mobeli can also help you if your tiles are smaller than 5.5" and do not have a non-porous surface? Two self-adhesive mobeli stainless steel plates are simply stuck on where you want them, provided that the base is solid. You don't have to drill or use any screws.



MOBELI®	Order no.
Plates of stainless steel/pair (5.7" x 5.7")	14614 14



MOBELI®

Mobility beyond limits. Portable grab bars.

Portable grab bar

Stay mobile. Overcome limitations.

The patented System Mobile from ROTH is a practical helper, which has proven itself for years in everyday life for disabled and elderly people.

The portable grab bars are quick and easy to fit and just as easy to remove. Both vacuum pads with a diameter of 4.6" can be pressed simply on a flat, non-porous surface, turn on the vacuum rocker – that's it.

All grab bars are made to a high standard and certified by the south German TÜV, have a strong adhesive force and are available in different lengths. Useful accessories such as, shower positioner, double hooks, tub shortener or tray – as well as numerous other helping "Hands" round up the System Mobeli sensibly.

MOBELI®	Grab length (inch)	Nominal load	Order no.
Dual grip	8.3	166 lbs	14002 25 S
Dual grip	13.8	154 lbs	14002 26 S
Dual grip	25.6	154 lbs	14002 27 S
Dual grip telescope	7.9-13.6	144 lbs	14002 21 S
Dual grip telescope	13.4-18.3	144 lbs	14002 22 S
Dual grip telescope	17.3-22.3	144 lbs	14002 23 S
Dual grip telescope	26.2-31.1	132 lbs	14002 24 S

MOBELI®



Vacuum pads

Made of special rubber. Also suitable for tiled surfaces in the bathroom which are slightly textured.

Telescopic handle

We developed this telescopic handle so that you can fix the grab bar regardless of the tile structure. The handles have by default a diameter of 1.4".

Quick locking

All Mobeli grab bar can be fixed without tools and easily dismantled for travelling. Like all metal parts, the quick lock is made of stainless steel.

Safety indicator

The safety indicator of the Mobeli grab bars is simply revolutionary. This is a world premiere: **the first grab bar which monitors the vacuum's strength and therefore the adhesive force.** The Mobeli safety rocker constantly monitors the vacuum. When the warning rocker is fully lowered, the adhesive force is at it's strongest and the grab bar is safe to use. If the warning rocker is raised showing the red warning rim, the Mobeli vacuum pads needs to be removed and reset. People with visibility impairments can touch the safety indicator and know by touching whether it's safe or not.

Important note: The Mobeli grab bars are available with or without the safety indicator. If you want the grab bar without the safety indicator, just leave the - S - off the order number.

Mobility Helpers

More mobility. On the go or at home.

High-performance version: The portable vacuum bar QuattroPlus can be used on the wall, in a corner or as a horizontal grip bar for training.

The integrated angle joint makes the QuattroPlus very flexible and its 4 vacuum pads increase the nominal load up to 276 lbs.

Diagonal joint adapters make it possible to fix the QuattroPlus in all directions and with a nominal load of up to 176 lbs.



MOBELI®	Grab length (inch)	Nominal load	Order no.
QuattroPlus	7.9-12.8	276 lbs	14002 41 S
QuattroPlus	12.6-17.5	276 lbs	14002 42 S
QuattroPlus	16.5-21.5	142 lbs	14002 43 S
QuattroPlus	25.4-30.3	142 lbs	14002 44 S

The field of application for all mobile telescopic grab bars can be enhanced considerably by clipping the cardan joint or angle joint adapter between the vacuum pad and the grip bar. Different angles can be achieved and a maximum of flexibility for positioning is made possible.

Mobeli is the original among the portable grab bars. Our product family is constantly being developed. ROTH, a German company rich in tradition, develops and manufactures the Mobeli grab bars in Germany and they are subject to German TÜV production monitoring.



MOBELI®	Nominal load	Order no.
Cardan joint adapter-pair	110 lbs	14005 55
Angle joint adapter-pair	depends on grab length	14005 58
Diagonal joint adapter-pair for QP	176 lbs	14615 57



Practical and comfortable. The tub shortener can be infinitely positioned and can prevent you from slipping down in the tub. It is also suitable as a tray in the bath.

Stabilizing handles. You want an aid that stabilizes your arm: Mobeli Stabi – in 2 versions. It can be fixed fast and easily to all non-porous surfaces such as varnished wood, plastic or metal.

The Stabi Horizontal stabilizes horizontally; Stabi Vario can be panned vertically, horizontally, or to the left or right at an angle of 45° as desired.



MOBELI®	Width of tub (inch)	Order no.
Bath tub shortener (2 suction pads)	from 20.1	14002 81 S
Bath tub shortener (3 suction pads)	20.1-27.6	14002 82 S

	Grip (inch)	Nominal load	Order no.
Stabi Horizontal	1.2	34 lbs	14002 95 S
Stabi Vario	1.2	34 lbs	14002 96 S

