

# Clarke Health Care Products

## Stroller Specifications/Comparison



	Corzo Basic 30	Corzo Basic 34	Corzo Basic 38	Corzo Xcountry 30	Corzo Xcountry 34	Corzo Xcountry 38	Corzo Xcountry 42	TDM Mini	TDM Standard	TDM Maxi	TDM Supermaxi
Product number	CNB30PYYY	CNB34PYYY	CNB38PYYY	CRX30CYYY	CRX34CYYY	CRX38CYYY	CRX42CYYY	T4WYCMTY	T4WYCYTY	T4WYCXTY	T4WYGSZY
Back height	20"	25.75"	33.5"	18" - 22"	22" - 25.75"	22" - 25.75"	26.5" - 30.5"	17" - 23.75"	20" - 25"	22" - 26"	22" - 26"
Seat depth	10"	15"	15"	10" - 13.25"	11.75" -	13"	15"	9.5" - 12"	9.5" - 12"	12.5" - 15.5"	12.5" - 16.75"
Seat width	11.5"	13.25"	15"	11.5"	13.25"	15"	16.5"	5" - 12"	9.5" / 11" / 12.5"	11" / 12.5" / 14.5"	13" / 15" / 16.25"
Seat to footrest	6.5"	10.5"	10.5"	6.5" - 14.25"	10.5" - 18.25"	10.5" - 18.25"	10.5" - 18.25"	5" - 12"	7.75" - 13.75"	9.25" - 15.5"	8.5" - 15"
Backrest angle	115°	115°	115°	115°	115°	115°	115°	0° - 180°	0° - 180°	0° - 180°	0° - 160°
Seat unit angle	25°	25°	25°	25°	25°	25°	25°	+/- 27°	+/- 27°	+/- 27°	+/- 27°
Total length	43.25"	44.5"	44.5"	44"	44.5"	44.5"	44.5"	33.5"	33.5"	36"	43"
Total width	21.5"	23"	24.75"	23"	24.5"	26"	38.5"	24"	24"	25.5"	27.5"
Total height w/ handle	36"	40"	40"	37.5"	38"	38"	30" - 40.75"	39" - 45"	39" - 45"	40" - 46"	42"
Weight	26.5 lbs	27.5 lbs	31 lbs	25 lbs	26.5 lbs	30 lbs	31 lbs	44 lbs	46.3 lbs	52 lbs	53 lbs
Capacity	100 lbs	120 lbs	165 lbs	100 lbs	120 lbs	165 lbs	165 lbs	77 lbs	77 lbs	110 lbs	165 lbs
Folded length	45"	48"	48"	45"	46"	46"	47"	36.5"	36.5"	37.5"	37.5"
Folded height	12"	14"	14"	14"	14.5"	13"	13"	21.5"	21.5"	24.5"	25.5"
Folded width	14"	14"	14"	15"	15"	15"	15"	24"	24"	25.25"	27.25"



Indigo
<b>LBU01</b>
21.5"
11"/13"/15"
11.5" - 14.5"
3.5" - 13.5"
0° - 25°
12° - 18° - 26°
na
na
44"
21.5 lbs
110 lbs
44"
12"
12.5"