

Jazz Specs

Jazz	Standard	Low
Basket capacity	10 pounds	10 pounds
Wheel	8"D x 1 ¼" w	8"D x 1 ¼" w
Capacity	287 lbs	287 lbs
Weight of walker	17 lbs	15 ¾ lbs
Overall width	24"	24"
Overall length	26"	26"
Width between wheels	17"	17"
Overall height	29 ½" – 37 ¾"	26" – 31 ½"
Seat height	23 ½"	20"
Seat	14 ½" L x 8" w	14 ½" L x 8"
Width between handles	15"	15"
Folded width	10"	10"