


# Clarke Health Care Products

## Walker Specifications Comparison



	Prod Number	Distance between rear wheels	Minimum handle height	Maximum handle height	Min & Max Tall handle height	Seat LxW	Seat height	Overall Width	Overall Length	Width between handles	Walker weight	Weight capacity
Legacy Standard	D12050	19"	33"	37"	35" - 39"	7" x 14"	24"	23"	26"	17"	15 lbs	350 lbs
Legacy Low	D12052	19"	31"	34"	34" - 37"	7" x 14"	22"	23"	26"	17"	15 lbs	350 lbs
Legacy Super Low	D12054	19"	28"	31"	31" - 35"	7" x 14"	19"	23"	23"	17"	15 lbs	350 lbs
Maxi+ Standard	D12120	22"	32"	38"	N/A	8.5" x 17"	25.5"	27.5"	30"	20.5"	21 lbs	440 lbs
Maxi+ Low	D12122	22"	28.5"	35"	N/A	8.5" x 17"	22"	27.5"	26"	20.5"	21 lbs	440 lbs
Symphony Standard	D12160	18.5"	30"	36"	N/A	7" x 14"	24"	27"	27.5"	17"	15 lbs	275 lbs
Symphony Low	D12161	18.5"	27.5"	33"	N/A	7" x 14"	21"	24"	24.5"	17"	15 lbs	275 lbs
Symphony Super Low	D12162	18.5"	25"	30"	N/A	7" x 14"	18"	23"	22"	17"	15 lbs	275 lbs
Alpha Advanced	D12220	20"	37.5"	47"	N/A	7" x 14"	24.5"	25.5"	27.5"	17"	28 lbs	275 lbs
Alpha Basic	D12230	20"	37.5"	47"	N/A	7" x 14"	24.5"	25.5"	27.5"	17"	28 lbs	275 lbs
Melody Standard	D12010	21"	30.5"	38"	N/A	7" x 14"	24"	24"	27"	17"	16.5 lbs	300 lbs
Melody Low	D12012	21"	27.5"	35"	N/A	7" x 14"	20.5"	24"	25"	16.5"	16.5 lbs	300 lbs
Melody Super Low	D12014	18"	25.5"	31"	N/A	7" x 14"	18.5"	21.5"	24.5"	14"	15.5 lbs	300 lbs
Soprano	D12030	20"	30.5"	37"	N/A	8.5" x 17"	23.5"	25.5"	33.5"	17"	21.5 lbs	355 lbs
Jazz 600	D12140	24"	29.25"	37.75"	N/A		23.5"		26"	15"	17 lbs	287 lbs
Jazz 520	D12142	24"	26"	31.5"	N/A		20"		26"	15"	15.75 lbs	287 lbs
Banjo	D1452441	17"	27.5"	38.25"	N/A	6.5" x 15"	23.75"	23.5"	27"	17"	19.75 lbs.	287 lbs.